

Trans Smart series

Steam sterilizers

TRANS
GETINGE
GETINGE GROUP



Trans – Getinge's new mid-range **value-for-money** offering

Getinge – the world's leading provider of solutions for infection control – is proud to introduce Trans, a complete new mid-range product line. They're based on the cornerstones of our high-end Getinge products, and deliver basic functionality at more affordable prices.

Trans products by Getinge are manufactured in our own facilities and comply with international standards, to ensure that you always get value for money.

Our long-term commitment to our customers and the full support from our global and well-trusted sales & service network make Trans by Getinge a truly reliable choice.



ISO 9001:2008



ISO 13485:2003



ISO 14001:2004



Trans Smart – steam sterilizers you can **rely on**



Trans by Getinge introduces Trans Smart Series Sterilizers that have been built to deliver top value-for-money performance.

Reliability

The Trans Smart series has been developed, tested and manufactured to relevant standards for steam sterilizers, and offer proven functionality and a robust design.

High capacity

Trans Smart series sterilizers give you surprisingly high loading capacity, with faster process times and greater throughput than other steam sterilizers in their price range.

Easy to use

Trans Smart series steam sterilizers are easy to use and easy to maintain. Their user-friendly HMI design includes a 7" Siemens PLC color touchscreen.

Simple and cost-effective

Trans Smart steam sterilizers are fast and versatile, and fulfill high demands on quality, efficiency and user-friendliness. They combine good capacity with cost-efficiency, while offering high operating reliability and easy maintenance.



New-generation touchscreens

Trans Smart features our new generation of user-friendly Siemens touchscreens with high resolution. The screens offer optimal ergonomics, with a clear overview from all angles in the control area.

The user interface and menu are intuitive, the text is remarkably distinct and the graphics are vivid.



Good quality 316 L

The acid-proof stainless steel chambers and steam jackets and enclosed in corrosion-proof aluminum to help assure long and trouble-free operation. CFC-free insulation is used throughout, including on appropriate piping.



Standard configuration

The vertically sliding door is fast and saves space, while offering proven, hands-free convenience.

The Trans Smart series comes with a built-in air compressor, electrical steam generator and wheels. It's also easy to install and move, making it highly versatile.





Easy access for service

All components are easily accessible for maintenance from both sides as standard (optional right- or left-side access only).

Short process time

The sterilization cycle in the new Trans Smart series is optimized for short process time and meets the European Norm EN285:2006 + A2:2009.

Clean steam

Provided that the integrated steam generators are supplied with treated water (the better the water quality, the better the steam). they will help prevent problems of soiling due to poor steam quality.



Versatile loading options

Every Trans Smart single-door steam sterilizer unit is supplied with a fixed-height loading trolley and a shelf rack. Each double-door unit is supplied with two fixed-height loading trolleys and a shelf rack.

Loading trolleys can be used for either baskets or direct placement of goods, and are equipped with a safety locking device. They are adapted to standardized basket and container systems.

Trans – Getinge's new mid-range product line – offers a complete, efficient, economical, ergonomical handling systems for baskets as well as containers. The systems may be used separately or in combination. And both systems are available with manual or loading and unloading. This means that your system can be tailored to your needs.



Manual loading/unloading of wire baskets



Manual loading/unloading of shelf rack



Manual loading of containers or various packs using a loading platform and loading module on a loading trolley. Also using optional loading module(s) to process containers



The shelf rack can also be loaded with wire baskets or other goods.



Loading side

Unloading side (double-door version)

MODEL	DOOR (S)	CODE	VOLUME (liters)	CHAMBER WxHxL (mm)	OUTSIDE DIMENSIONS WxHxL (mm)	30x30x60 CONTAINER QUALITY (STU)	600x400x200 CONTAINER QUALITY (ISO)
Smart	Single	TS130EC1NX	130	500x500x540	900x1800x1100	1/2	0
Smart	Single	TS160EC1NX	160	350x670x700	770x1950x1300	2	2
Smart	Single	TS180EC1NX	180	500x500x724	900x1800x1300	1	1
Smart	Single	TS250EC1NX	250	500x500x1000	900x1800x1600	1	1
Smart	Single	TS330EC1NX	330	670x670x742	1050x1950x1300	4	2
Smart	Single	TS420EC1NX	420	670x670x943	1050x1950x1600	6	4
Smart	Single	TS555EC1NX	555	670x670x1245	1050x1950x1800	8	6
Smart	Single	TS670EC1NX	670	670x670x1496	1050x1950x2000	8	6
Smart	Single	TS740EC1NX	740	670x670x1655	1050x1950x2250	10	8
Smart	Single	TS825EC1NX	825	670x670x1842	1050x1950x2400	12	8
Smart	Single	TS900EC1NX	900	630x900x1601	1500x1950x2200	10	6
Smart	Single	TS1000EC1NX	1000	630x900x1770	2200x1750x2000	10	6
Smart S	Double	TS160EC2NX	160	350x670x700	770x1950x1200	2	2
Smart S	Double	TS250EC2NX	250	500x500x1000	1090x1800x1500	1	1
Smart S	Double	TS330EC2NX	330	670x670x742	1250x1950x1250	4	2
Smart S	Double	TS420EC2NX	420	670x670x943	1250x1950x1440	6	4
Smart S	Double	TS555EC2NX	555	670x670x1245	1250x1950x1745	8	6
Smart S	Double	TS670EC2NX	670	670x670x1496	1250x1950x2000	8	6
Smart S	Double	TS740EC2NX	740	670x670x1655	1250x1950x2160	10	6
Smart S	Double	TS825EC2NX	825	670x670x1842	1250x1950x2345	12	8
Smart S	Double	TS900EC2NX	900	630x900x1601	1500x1950x2100	10	6
Smart S	Double	TS1000EC2NX	1000	630x900x1770	2200x1750x2000	10	6

* Basket size (SPRi/ISO): 585x395x195 / 600x400x200 mm (LxWxH)

** STU size: 600x300x300 mm (LxWxH)



Trans by Getinge provides value for money sterilizers that are compliant to international standards. Our service network gives you reliable support you can depend on for your daily operations. Trans by Getinge is a brand of the Getinge Group.

www.transmedikal.com

GETINGE
GETINGE GROUP

Legal Manufacturer
TRANS MEDİKAL ALETLER SAN. VE TIC. A.Ş
Ostim(OSB) Mah. 1219. Sok. No:1
Yenimahalle- Ankara- Türkiye
Tel: + 90 (312) 385 77 20-36
Fax: + 90 (312) 385 25 50
info@transmedikal.com
www.transmedikal.com

GETINGE GROUP

Getinge Group is a leading global provider of products and systems that contribute to quality enhancement and cost efficiency within healthcare and life sciences. We operate under the three brands of ArjoHuntleigh, Getinge and Maquet. ArjoHuntleigh focuses on patient mobility and wound management solutions. Getinge provides solutions for infection control within healthcare and contamination prevention within life sciences. Maquet specializes in solutions, therapies and products for surgical interventions, interventional cardiology and intensive care.
